

Product Code . DS-C-11341

Test Tube Holder



Description

Catalogue No. : "IDAC1579/25"

DESCRIPTION :-

- Material: wooden
- A test tube holder is used to hold test tubes.
- It is used for holding a test tube in place when the tube is hot or cannot be touched.
- With full accuracy.
- Used to hold the test tube when the test tube is being heated.

Didac Scientific Pvt. Ltd,
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom
Direct Contact Details ✉ sales@didacscientific.co.uk
💻 www.didacscientific.co.uk/