



Product Code . DS-C-11345

## Thistle Funnel

### Description

---

Catalogue No. : "IDAC5299/29"

#### DESCRIPTION :-

- Material: Borosilicate Glass
- Length : 300 mm
- Used in gas preparation, osmosis etc.
- Chemical resistance
- Color : Clear
- Heat resistance
- Color: transparent

---

Didac Scientific Pvt. Ltd,  
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom  
Direct Contact Details ✉ [sales@didacscientific.co.uk](mailto:sales@didacscientific.co.uk)  
💻 [www.didacscientific.co.uk/](http://www.didacscientific.co.uk/)