



Product Code . DS-C-11346

## Red Litmus Paper strip

### Description

---

Catalogue No. : "IDAC6501/30"

#### DESCRIPTION :-

- Red – color changes from red to blue when dissolved in base solution.
- Litmus indicator solution turns red in acidic solutions and blue in alkaline solutions - and it turns purple in neutral solutions.
- Litmus paper is usually more reliable, and comes as red litmus paper and blue litmus paper.
- The table shows the color changes it can make.
- 100 leaves as per pack

---

Didac Scientific Pvt. Ltd,  
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom  
Direct Contact Details ✉ [sales@didacscientific.co.uk](mailto:sales@didacscientific.co.uk)  
💻 [www.didacscientific.co.uk/](http://www.didacscientific.co.uk/)