



Product Code . DS-C-11347

Blue Litmus Paper Strip

Description

Catalogue No. : "IDAC7205/31"

DESCRIPTION :-

- Blue – color changes from Blue to red when dissolved in base solution.
- The blue litmus will turn red or pink when it comes into contact with an acidic solution below pH 4.5.
- Between 4.5 and 8.3 the paper will be purple.
- Above 8.3 it will stay blue.
- Supplied in plastic vial, about 100 paper strips/vial.
- 100 leaves as per pack

Didac Scientific Pvt. Ltd,
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom
Direct Contact Details ✉ sales@didacscientific.co.uk
🌐 www.didacscientific.co.uk/