



Product Code . DS-C-11350

Boss Head S type

Description

Catalogue No. : "IDAC5466/32.2"

DESCRIPTION :-

- Double ended
- Material: Aluminum
- For holding square/round rod perpendicular to each other.
- Clamp opening 16mm.
- 't' thumb screws made of chromed iron.

Didac Scientific Pvt. Ltd,
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom
Direct Contact Details ✉ sales@didacscientific.co.uk
💻 www.didacscientific.co.uk/