



Product Code . DS-C-11351

Test Tube Rack

Description

Catalogue No. : "IDAC4729/33"

DESCRIPTION :-

- Test tube stand plastic 6 holes 6 pegs
- This autoclavable stand is perhaps the most commonly used test tube stand in various laboratories. For better visibility, this 6place rack is end supported by 2 columns.
- The top plate has 2 holes for holding 25 mm dia test tubes, four holes to hold 16 mm dia tubes where as the base has corresponding hemispherical wells.
- It also has six vertical pins for drying tubes. Another option available in this section is a test tube stand in which all the holes in the top
- Plates can hold 25 mm dia tubes.

Didac Scientific Pvt. Ltd,
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom
Direct Contact Details ✉ sales@didacscientific.co.uk
💻 www.didacscientific.co.uk/