



Product Code . DS-C-11364

Voltmeter

Description

Catalogue No. : "IDAC2999/04"

DESCRIPTION :-

- DC range: 0-10V
- A voltmeter, also known as a voltage meter, is an instrument used for measuring the potential difference, or voltage, between two points in an electrical or electronic circuit.
- It is generally used in all types of Laboratories such as School, college, university, Research etc.
- Desk Stand: Top Terminal ABS plastic moulded, Meters moving coil dc.

Didac Scientific Pvt. Ltd,
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom
Direct Contact Details ✉ sales@didacscientific.co.uk
💻 www.didacscientific.co.uk/