



**Product Code . DS-C-11365**

## **Electronic Balance**

### **Description**



---

Catalogue No. : "IDAC3102/05"

#### **DESCRIPTION :-**

- Capacity: 620gm
- Die-cast aluminum housing
- Variety of measurement units
- Piece-counting
- Percent weighing
- Density determination
- Accuracy: 0.01g, Readability - 0.01g With Wind Shield.
- Zero tracking, Low Battery Indication.
- Rechargeable Battery.

- 
- Adopter
  - Power Supply 220/230V, 50/60Hz
- 

Didac Scientific Pvt. Ltd,  
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom  
Direct Contact Details  [sales@didacscientific.co.uk](mailto:sales@didacscientific.co.uk)  
 [www.didacscientific.co.uk/](http://www.didacscientific.co.uk/)