

Product Code . DS-C-11371

Conical Flask



Description

Catalogue No. : "IDAC9311/11"

DESCRIPTION :-

- Capacity : 250ml
- Made of Polypropylene, these cone shaped flasks are rigid, translucent & autoclavable.
- The screw cap keeps the contents of the flask contamination free and makes the flask leakproof.
- Threads be disengaged before autoclaving.

Didac Scientific Pvt. Ltd,
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom
Direct Contact Details ✉ sales@didacscientific.co.uk

