

Product Code . DS-C-11374

## Conical Flask



### Description

---

Catalogue No. : "IDAC3922/14"

DESCRIPTION :-

- Capacity: 250ml
- Material: Polypropylene
- Made of Polypropylene, these cone shaped flasks are rigid, translucent & autoclavable.
- The screw cap keeps the contents of the flask contamination free and makes the flask leakproof.
- Threads be disengaged before autoclaving.

---

Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom

Direct Contact Details  [sales@didacscientific.co.uk](mailto:sales@didacscientific.co.uk)

 [www.didacscientific.co.uk/](http://www.didacscientific.co.uk/)