

Product Code . DS-C-11380

Burette Stand



Description



Catalogue No. : "IDAC9113/20"

FEATURES:

- Best quality raw material
- Long life
- Most appreciated

DESCRIPTION :-

- Size : 225 X 150mm
- With side hole
- These non-corrosive, moulded Retort Stands provide excellent replacement for the metal ones.
- Even the metallic rod has a polypropylene covering making it non-corrosive. The rod measuring
- 12.5 mm x 750 mm can be screwed into the heavy base giving way to a perfectly stable Retort Stand.

Didac Scientific Pvt. Ltd,
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom
Direct Contact Details  sales@didacscientific.co.uk
 www.didacscientific.co.uk/