

Product Code . DS-C-11381

Fisher Clamp single



Description

Catalogue No. : "IDAC1105/20.1"

DESCRIPTION :-

- FISHER CLAMP SINGLE
- These Clamp are moulded in Polypropylene making it completely corrosion proof.
- The spring loaded jack in the middle provides a vice like grip to the Burette and prevents it from slipping.
- Two tapered wedges are provided to ensure firm grip on rods having dia
- $\frac{1}{2}$ " or $\frac{3}{8}$ ".

Didac Scientific Pvt. Ltd,
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom
Direct Contact Details ✉ sales@didacscientific.co.uk
🌐 www.didacscientific.co.uk/