

Product Code . DS-C-11382

Burette





Description

Catalogue No. : "IDAC2117/20.2"

DESCRIPTION :-

- Capacity : 25ml
- Material: Borosilicate Glass
- With Stopcock
- Heat Resistant
- Chemical Resistant
- Graduated
- This Burette has a transparent acrylic body and a leakproof stopcock fitted with a self lubricating PTFE plug.
- Each Burette is individually calibrated to give maximum accuracy.
- Acrylic body is unaffected by dilute mineral acids & bases except Hydrofloric Acid & Ammonium Hydroxide.

Direct Contact Details  sales@didacscientific.co.uk
 www.didacscientific.co.uk/