

Product Code . DS-C-11388

Wash bottle





Description

Catalogue No. : "IDAC8529/26"

DESCRIPTION :-

- Material : Polythene
 - Capacity : 250ml
 - Wash Bottles are made of Low Density Polyethylene which gives these bottles a translucent & unbreakable character.
 - Solutions can be poured easily by squeezing the bottle lightly.
 - The cap is fitted with a flexible polyethylene delivery tube which can be aimed wherever required.
-

Didac Scientific Pvt. Ltd,
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom
Direct Contact Details  sales@didacscientific.co.uk
 www.didacscientific.co.uk/