



Product Code . DS-C-11398

Buchner Funnel

Description

Catalogue No. : "IDAC9865/36"

DESCRIPTION :-

- Diameter : 70mm
- Moulded in Polypropylene, these two-piece funnels are light weight, easy to clean & autoclavable.
- The top portion can be snap fitted with the lower one forming a vacuum seal which does not require greasing.
- The base of the top portion has a perforated filter plate.?

Didac Scientific Pvt. Ltd,
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom
Direct Contact Details ✉ sales@didacscientific.co.uk
💻 www.didacscientific.co.uk/