



Product Code . DS-C-11401

Beaker

Description

Catalogue No. : "IDAC3479/39"

DESCRIPTION :-

- Material : Polypropylene
- Capacity : 50ml
- Color : White
- Graduated
- Beakers, moulded in Polypropylene, have excellent clarity and very good chemical resistance.
- These autoclavable Beakers have prominently raised graduations on the outer surface.
- Pouring a solution becomes very convenient because of its tapered spout.

Didac Scientific Pvt. Ltd,
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom
Direct Contact Details ✉ sales@didacscientific.co.uk
🌐 www.didacscientific.co.uk/