



Product Code . DS-C-11402

## Calorimeter

### Description



---

Catalogue No. : "IDAC4825/40"

DESCRIPTION :-

- Inner vessel: Made of Copper
  - Size: 5.1 cm x 7.5 cm (dia. x height)
  - Outer vessel: Made of Copper (electroplated) Size: 7 cm x 9.8 cm (dia. x height)
  - Material of cover Lid: Bakelite
  - Material: Copper
  - Good quality, polished finish vessel.
  - With wooden/bake lite lid fitted.
  - With rubber stopper and metal stirrer.
-

---

Didac Scientific Pvt. Ltd,  
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom  
Direct Contact Details  [sales@didacscientific.co.uk](mailto:sales@didacscientific.co.uk)  
 [www.didacscientific.co.uk/](http://www.didacscientific.co.uk/)