



Product Code . DS-C-11403

Stirrer

Description

Catalogue No. : "IDAC5969/41"

DESCRIPTION :-

- Material : Polypropylene
- Color : White
- Chemical resistant
- Made of polypropylene, these Stirrers are very useful for stirring different solutions in the laboratory.
- These Stirrers prove to be tough and long lasting.
- Dimension : Diameter 10mm x Height 250mm

Didac Scientific Pvt. Ltd,
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom
Direct Contact Details ✉ sales@didacscientific.co.uk
🌐 www.didacscientific.co.uk/