

Product Code . DS-C-11404

Thermometer



Description

Catalogue No. : "IDAC6196/42"

DESCRIPTION :-

- Type : Alcohol
- It is the general purpose stirring thermometer
- Length 300mm
- TEMP RANGE -10 TO 110°C
- SUBDIVISION DEGREES 1.

Didac Scientific Pvt. Ltd,
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom
Direct Contact Details ✉ sales@didacscientific.co.uk
💻 www.didacscientific.co.uk/