

Product Code . DS-C-11407

## Separating funnel



### Description



---

Catalogue No. : "IDAC0205/01"

DESCRIPTION :-

- SEPARATING FUNNEL 250ml
- These Separatory Funnel is made of Polypropylene, giving it excellent contact clarity and strength.

- 
- This pear shaped, autoclavable separating funnel has a leakproof threaded screw cap.
  - The Polypropylene stopcock, fitted with a self lubricating PTFE plug, is push fitted to the stem of the separatory funnel
- 

Didac Scientific Pvt. Ltd,  
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom  
Direct Contact Details  [sales@didacscientific.co.uk](mailto:sales@didacscientific.co.uk)  
 [www.didacscientific.co.uk/](http://www.didacscientific.co.uk/)