

Product Code . DS-C-11410

Stand with Clamp



Description

Catalogue No. : "IDAC4826/04"



DESCRIPTION :-

- Retort Stand Base (without Rod) 200x125mm
- Retort Stand base made of mild steel sheet Powder coated finish.
- Supplied without rod.

-
- A Retort Stand is a piece of scientific equipment to which clamps can be attached to hold other pieces of equipment and glassware for instance, burettes, test tubes and flasks.
 - Made of cast iron powder coated mild steel rod chrome plated.

Features:

- Best quality raw material
- Long life
- Most appreciated

Didac Scientific Pvt. Ltd,
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom
Direct Contact Details  sales@didacscientific.co.uk
 www.didacscientific.co.uk/