



Product Code . DS-C-11412

Boss Head S type

Description

Catalogue No. : "IDAC67234/4.2"

DESCRIPTION :-

- Boss Head S Type
- Double ended, made of Aluminum for holding square/round rod perpendicular to each other.
- Clamp opening 16mm.
- High Efficiency
- Anti-corrosive
- 'T' thumb screws made of chromed Iron.

Didac Scientific Pvt. Ltd,
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom
Direct Contact Details ✉ sales@didacscientific.co.uk
💻 www.didacscientific.co.uk/