



Product Code . DS-C-11413

Red Litmus Paper

Description

Catalogue No. : "IDAC7239/05"

DESCRIPTION :-

- Litmus indicator solution turns red in acidic solutions and blue in alkaline solutions - and it turns purple in neutral solutions.
- Litmus paper is usually more reliable, and comes as red litmus paper and blue litmus paper.
- The table shows the color changes it can make.
- Red – color changes from red to blue when dissolved in base solution.

Didac Scientific Pvt. Ltd,
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom
Direct Contact Details ✉ sales@didacscientific.co.uk
💻 www.didacscientific.co.uk/