

Product Code . DS-C-11414

## Blue Litmus Paper



### Description

---

Catalogue No. : "IDAC8241/06"

#### DESCRIPTION :-

- The blue litmus will turn red or pink when it comes into contact with an acidic solution below pH 4.5.
- Between 4.5 and 8.3 the paper will be purple.
- Above 8.3 it will stay blue.
- Supplied in plastic vial, about 100 paper strips/vial.

---

Didac Scientific Pvt. Ltd,  
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom  
Direct Contact Details ✉ [sales@didacscientific.co.uk](mailto:sales@didacscientific.co.uk)  
🌐 [www.didacscientific.co.uk/](http://www.didacscientific.co.uk/)