



**Product Code . DS-C-11418**

## **Tripod Stand**

### **Description**



---

Catalogue No. : "IDAC3267/10"

#### **DESCRIPTION :-**

- Height x side of triangular 15x12.5cm
- Material: Cast iron top and mild steel legs painted.
- Zinc coated legs are bent outwards for stability.
- A tripod is a portable three-legged frame, used as a platform for supporting the weight and maintaining the stability of some other object.
- A tripod provides stability against downward forces and horizontal forces and movements about horizontal axes.

- 
- A laboratory tripod is a piece of three-legged equipment commonly used to conduct experiments in laboratories.
  - It is used as a platform to hold and support glassware, such as beakers and flasks, during experiments
- 

Didac Scientific Pvt. Ltd,  
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom  
Direct Contact Details  [sales@didacscientific.co.uk](mailto:sales@didacscientific.co.uk)  
 [www.didacscientific.co.uk/](http://www.didacscientific.co.uk/)