



Product Code . DS-C-11420

Bunsen Burner

Description

Catalogue No. : "IDAC8280/12"

DESCRIPTION :-

- BURNER BUNSEN WITHOUT STOPCOCK
- Nickel plated burner tube of 12mm diameter with air regulator and a riffled connector.
- Mounted on a painted cast iron base.
- Only the air supply port can be adjusted.
- Without stopcock

Didac Scientific Pvt. Ltd,
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom
Direct Contact Details ✉ sales@didacscientific.co.uk
🌐 www.didacscientific.co.uk/