

Product Code . DS-C-11425

Magnet



Description

Catalogue No. : "IDAC2302/20"

DESCRIPTION :-

- Size : 100mm
- Material : Alnico
- This pair of plastic cased bar magnets is perfect for teaching magnetic attraction and repulsion.
- The magnets are durable and colored for ease of identifying north and south

Didac Scientific Pvt. Ltd,
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom
Direct Contact Details ✉ sales@didacscientific.co.uk
💻 www.didacscientific.co.uk/