

Product Code . DS-C-11430

Big Hard Glass Test Tube



Description

Catalogue No. : "IDAC1310/26"

DESCRIPTION :-

- TEST TUBE BOROSILICATE GLASS WITH RIM 32X200mm
- Material : Borosilicate Glass
- With Rim
- Heat resistant
- Chemical resistant
- Color : Clear
- Size : 32x200mm

Didac Scientific Pvt. Ltd,
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom
Direct Contact Details  sales@didacscientific.co.uk
 www.didacscientific.co.uk/