

Product Code . DS-C-11432

Funnel



Description

Catalogue No. : "IDAC2315/28"

DESCRIPTION :-

- Material : Polypropylene
- Designed specifically for analytical chemistry, the body of the funnels is 600.
- Moulded in PP, these funnels are provided with internal & external ribs to prevent air blockage.
- Diameter 75mm. Inside Dia: 75mm
- Stem Length: 77mm

Didac Scientific Pvt. Ltd,
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom
Direct Contact Details ✉ sales@didacscientific.co.uk
🌐 www.didacscientific.co.uk/