

Product Code . DS-C-11433

Thermometer



Description

Catalogue No. : "IDAC3319/29"

DESCRIPTION :-

- Type : Alcohol
- It is the general-purpose stirring thermometer
- Length : 76mm
- Length : 300mm
- Temp Range : -10 TO 110°C
- Subdivision Degrees 1.

Didac Scientific Pvt. Ltd,
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom
Direct Contact Details ✉ sales@didacscientific.co.uk
💻 www.didacscientific.co.uk/