



Product Code . DS-C-11436

Rubber Stopper

Description

Catalogue No. : "IDAC8522/32"

DESCRIPTION :-

- Two hole no.2
- Top diameter : 16.5mm
- Bottom diameter : 13.0mm
- Height 20.0mm
- Material : Rubber
- Chemical resistant

Didac Scientific Pvt. Ltd,
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom
Direct Contact Details ✉ sales@didacscientific.co.uk
💻 www.didacscientific.co.uk/