



Product Code . DS-C-11440

Rubber Stopper

Description

Catalogue No. : "IDAC9613/32.4"

DESCRIPTION :-

- Two hole no.6
- Top diameter : 24.5mm
- Bottom diameter : 20.0mm
- Height : 24.0mm.
- Material : Rubber
- Chemical resistant

Didac Scientific Pvt. Ltd,
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom
Direct Contact Details ✉ sales@didacscientific.co.uk
💻 www.didacscientific.co.uk/