

Product Code . DS-C-11443

Conical Flask



Description

Catalogue No. : "IDAC4519/36"

DESCRIPTION :-

- Capacity: 250ml
- Material: Polypropylene
- Made of Polypropylene, these cone shaped flasks are rigid, translucent & autoclavable.
- The screw cap keeps the contents of the flask contamination free and makes the flask leakproof.
- Threads be disengaged before autoclaving.

Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom

Direct Contact Details  sales@didacscientific.co.uk

 www.didacscientific.co.uk/