

Product Code . DS-C-11445

Measuring cylinders





Description

Catalogue No. : "IDAC4893/38"

DESCRIPTION :-

- Material : Polypropylene
- The Measuring Cylinders, moulded in Polypropylene & TPX, are single piece cylinders.
- These autoclavable cylinders have raised prominent graduations for easy reading and have excellent contact clarity.
- The broad base provides good stability.
- Capacity: 250ml
- Least count: 2ml?

Didac Scientific Pvt. Ltd,
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom
Direct Contact Details  sales@didacscientific.co.uk
 www.didacscientific.co.uk/