



Product Code . DS-C-11450

Rubber Stopper

Description

Catalogue No. : "IDAC9335/40.3"

DESCRIPTION :-

- Rubber stopper hole no.5
- Top diameter : 22.5mm
- Bottom diameter : 17.5mm
- Height : 23.0mm.
- Material : Rubber
- Chemical resistant

Didac Scientific Pvt. Ltd,
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom
Direct Contact Details ✉ sales@didacscientific.co.uk
💻 www.didacscientific.co.uk/