



Product Code . DS-C-11453

Atomic Model Set

Description

Catalogue No. : "IDAC1147/41"

DESCRIPTION :-

- Containing senior 370 balls
- These sets consist of moulded balls of different colours and sizes along with connecting lugs of different sizes.
- These sets are very helpful in modelling different organic Inorganic compounds.
- This set consists of 150 lugs 370 balls of prescribed colours and sizes packed in a moulded box having separate space for every particular colour size of ball.

Didac Scientific Pvt. Ltd,
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom
Direct Contact Details ✉ sales@didacscientific.co.uk
💻 www.didacscientific.co.uk/