

Product Code . DS-C-11454

## Four Finger Clamp





### Description

---

Catalogue No. : "IDAC851/42"

DESCRIPTION :-

- 
- Made of Aluminum-alloy die casting.
  - This Laboratory clamp has PVC coated jaws capable of holding upto 90 mm article.
  - It is provided with 8mm diameter & 6" long Stain less Steel rod.
  - It is coated with LEAD-FREE powder coating & fitted with brass screw.
- 

Didac Scientific Pvt. Ltd,  
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom  
Direct Contact Details  [sales@didacscientific.co.uk](mailto:sales@didacscientific.co.uk)  
 [www.didacscientific.co.uk/](http://www.didacscientific.co.uk/)