



Product Code . DS-EE-11477

Voltmeter

Description

AC Voltmeter

Principle : Moving coil type

Accuracy : Class 0.5

Scale Length : 135mm (5-3/8inch)

Weight : 1.7kg

Dimension : 195(Length) x 169(Width) x 86(Depth) mm

Case Material: Bakelite with glass window

4 Ranges:75/150/300/600V

Features :-

High efficiency and convenient use

Compact size

Easy operation and stable work

Binding post or safety socket connection

DC Voltmeter

Principle : Moving coil type

Accuracy : Class 0.5

Scale Length : 135mm (5-3/8inch)

Weight : 1.7kg

Dimension : 195(Length) x 169(Width) x86(Depth) mm

Case Material: Bakelite with glass window

4 Ranges: 1v-1000v

Features :-

High efficiency and convenient use

Compact size

Easy operation and stable work

Binding post or safety socket connection

Didac Scientific Pvt. Ltd,
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom
Direct Contact Details  sales@didacscientific.co.uk
 www.didacscientific.co.uk/