

Product Code . DS-EE-11489

AC & DC Generation Model



Description

Technical Features:-

- Learn about current and see a visual demonstration
- The elegant design produces AC & DC current at the same time
- Allows students to easily compare and contrast an AC generator and a DC generator
- Activity guide included
- AC Motor/Generator
- Motor: AC 220V, 480W, 1500rpm, 12V DC EXV
- Generator: AC 220V, 320W, 1500rpm
- Prime Mover: 24V DC, 1200mA, 1500rpm
- DC Motor/Generator
- Motor: 24V DC, 1200mA, 1500rpm
- Generator: 24V DC, 750mA, 1500rpm
- Prime Mover: AC 220V, 480W, 1500rpm, 12V DC EXV

