



Product Code . DS-EE-11496

Galvanometer

Description

Technical Features:-

- These low-cost devices are designed for use as desk-mounted / table-top units and have rotary winding units, mirror scales and safety sockets in 4mm technology.
- These analogue devices do not require any voltage supply, which makes them durable and cost-saving to use.
- On the large analogue mirror scale, the respective measured values can be conveniently read by the user or student.
- Compact size
- Easy operation and stable work
- Binding post or safety socket connection
- Type: Galvanometer
- Range : $\pm 35 \mu\text{A}$