



Product Code . DS-EE-11522

## Prony Breck

### Description

---

Catalogue No. : "IDAC127/012"

#### DESCRIPTION :-

- Material: Stainless steel
- Complete working model with set of weights and technical manual
- Rust Proof parts ensure long life
- Long Grips for easy carrying around
- Thickness: 8mm , mounted on ivory Bakelite board.
- Steel grip: 150mm
- Fly wheel: 100mm,
- Torque arm: 500mm, with weight set & spring balance.

---

Didac Scientific Pvt. Ltd,  
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom  
Direct Contact Details ✉ [sales@didacscientific.co.uk](mailto:sales@didacscientific.co.uk)  
💻 [www.didacscientific.co.uk/](http://www.didacscientific.co.uk/)