



Product Code . DS-EE-11542

Outside Micrometer

Description

Catalogue No. : “IDAC205/032”

DESCRIPTION :-

- Outside micrometer has a vernier scale in metric units for taking precise outside diameter (OD) measurements, and a ratchet-stop mechanism for uniform pressure during adjustment.
- Measuring range- 25-50 mm.
- Additional scale-0.01 mm
- Operating temperature-5 Degree Celsius - 40 Degree Celsius
- Influence of humidity- Not significant below 80%
- Material- Aluminum Carbon Steel.

Didac Scientific Pvt. Ltd,
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom
Direct Contact Details ✉ sales@didacscientific.co.uk
💻 www.didacscientific.co.uk/