



Product Code . DS-EE-11567

Working Table

Description

Catalogue No. : "IDAC315/057"

DESCRIPTION :-

- Table Top Material : Garjan Wood mm 37.5 mm
- Table leg/ frame materials: Mild Steel 50mm x 50mm x 5mm angle bar
- Table Height : 760mm
- Table length Size : 2400 mm
- Table Width Size : 1200 mm

Didac Scientific Pvt. Ltd,
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom
Direct Contact Details ✉ sales@didacscientific.co.uk
💻 www.didacscientific.co.uk/