

Product Code . DS-EE-11576

Volt Meter



Description

Catalogue No. : "IDAC113/065"

DESCRIPTION :-

- Range: A/C 250V
- Deflection Angle: 90°
- Scale Length : 65mm
- DC Voltmeter:
 - Range: D/C 3/30/300V
 - Deflection Angle: 90°
 - Scale Length: 65mm
- Terminal: 4mm socket captive head

Didac Scientific Pvt. Ltd,
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom
Direct Contact Details ✉ sales@didacscientific.co.uk
💻 www.didacscientific.co.uk/