



Product Code . DS-LE-11586

Galvanometers

Description

Catalogue No. : "IDAC109/03"

DESCRIPTION :-

Centre Zero Universal Meter Galvanometer

The centre zero universal meter (Universal Galvanometer) is electrically identical to the universal meter.

1000 ohm movement, moving coil, diode protection against overload up to 12V D.C.

Didac Scientific Pvt. Ltd,
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom
Direct Contact Details ✉ sales@didacscientific.co.uk
💻 www.didacscientific.co.uk/