



**Product Code . DS-LE-11587**

## **Triple Beam Balance**

### **Description**

---

**Catalogue No. : “IDAC111/04”**

#### **DESCRIPTION :-**

Triple Beam Balance are specially designed to handle weighing solids and powder.

Three notched and tiered beams with center reading.

Indicating sliding masses giving total weighing capacity of 610 gm.

The capacity can be increased up to 2610 gm. by using supplementary masses.

Spring loaded zero adjust compensator and magnetic damping are provided for fast setting and reading..

Capacity : 610 gm.

Readability : 0.1 gm.

Front Beam Calibration : 10 x 0.1 gm.

Second Beam Calibration : 500 x 100 gm.

Rear Beam Calibration : 100 x 10 gm.

---

Didac Scientific Pvt. Ltd,  
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom  
Direct Contact Details ✉ [sales@didacscientific.co.uk](mailto:sales@didacscientific.co.uk)  
💻 [www.didacscientific.co.uk/](http://www.didacscientific.co.uk/)