

Product Code . DS-LE-11588

Electronic Balance



Description

Catalogue No. : "IDAC113/05"

DESCRIPTION :-

Readability : 0.001g

Die-cast aluminum housing

Variety of measurement units

Piece-counting

Percent weighing

Density determination

Didac Scientific Pvt. Ltd,
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom
Direct Contact Details ✉ sales@didacscientific.co.uk
💻 www.didacscientific.co.uk/