



Product Code . DS-LE-11592

Bee Hive Shelves

Description

Catalogue No. : "IDAC121/09"

DESCRIPTION :-

Used to support a receiving jar or tube while a gas is being collected over water with a pneumatic trough.

Material: Porcelain

Used in Laboratory

Color : White

Didac Scientific Pvt. Ltd,
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom
Direct Contact Details ✉ sales@didacscientific.co.uk
💻 www.didacscientific.co.uk/