



Product Code . DS-LE-11596

Vernier Caliper

Description

Catalogue No. : “IDAC129/013”

DESCRIPTION :-

A Vernier Caliper is a visual aid to take an accurate measurement reading between two graduation markings on a linear scale by using mechanical interpolation; thereby increasing resolution and reducing measurement uncertainty by using Vernier acuity to reduce human estimation error.

Plated steel with dual scale, 0-15cm x 0.1mm and 6 x 1/1283.

Provided with thumb movement for easy motion of jaws.

Didac Scientific Pvt. Ltd,
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom
Direct Contact Details ✉ sales@didacscientific.co.uk
💻 www.didacscientific.co.uk/