

Product Code . DS-LE-11599

Retort Stands (Clamp Stands Metal) Metal



Description

Catalogue No. : “IDAC140/016”

DESCRIPTION :-

A Retort Stand is a piece of scientific equipment to which clamps can be attached to hold other pieces of equipment and glassware for instance, burettes, test tubes and flasks.



Made of cast iron powder coated mild steel rod chrome plated.

Features:

Best quality raw material

Long life

Most appreciated

Didac Scientific Pvt. Ltd,
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom
Direct Contact Details  sales@didacscientific.co.uk
 www.didacscientific.co.uk/